



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® round bottom flask, wide neck, with beaded rim, 50 ml

€82.30
plus VAT &
Shipping

Product Images



Description

DURAN® round-bottomed flasks, wide-necked, not according to DIN EN ISO, 50 ml Bellied vessels are called flasks in the laboratory sector and, in addition to beakers, are other universally applicable vessels. They are used, for example, for distillation or for heating liquids. Flasks are available in a wide variety of designs: narrow- and wide-necked, with or without ground joint, as round-bottomed or upright flasks. According to the application, the required vessel with its specific components is selected. Special features and advantages: The wide neck enables easy filling and cleaning. The geometry enables very uniform heating. Burned-in and therefore permanently durable white printed image. Application: use in heating mantles, heating, mixing of several components, sizes 50, 1000, 2000, 3000, 5000, 6000, 10000, 12000 and 20000 ml not according to DIN EN ISO.

Additional Information

No.	SD-0746
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112174117
EAN	4058072181826
Manufacturer EAN	4032051013513
Gross weight	0.05
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	0.04000 kg
Material	Borosilicate glass 3.3
Volume	50 ml
sterile	No
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packaging volume	0.00331800 cc
Customs tariff number	70172000
Number of side necks	1
Neck size	far

