



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz three-neck round bottom flask, side necks
slanted, 100 ml, MH NS 14/23, SH NS 14/23**

**€41.76
plus VAT &
Shipping**

Product Images

Description

Reaction flasks according to DIN 12392 or 12394 and in addition to DIN, with two slanted (angle: 20°) side necks, with standard ground joint and with labeling field, material: DURAN® (borosilicate glass 3.3)

Additional Information

No.	LZ-0667
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2667
EAN	4058072304072
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Pistons properties	with side necks with standard grinding
Number of side necks	3
Alignment side necks	Oblique 20°
DM Piston	64 mm
Standard ground piston	14/23
Standard grinding for side necks	14/23
Type piston	Three-neck round-bottom flask

