



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz three-neck round bottom flask, side necks slanted, 20,000 ml, MH NS 45/40, SH NS 29/32

€947.66
plus VAT &
Shipping

Product Images



Description

Reaction flasks according to DIN 12392 or 12394 and in addition to DIN, with two slanted (angle: 20°) side necks, with standard ground joint and with labeling field, material: DURAN® (borosilicate glass 3.3)

Additional Information

No.	LZ-0707
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2707
EAN	4058072304478
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	20000 ml
Pistons properties	with side necks with standard grinding
Number of side necks	3
Alignment side necks	Oblique 20°
DM Piston	345 mm
Standard ground piston	45/40
Standard grinding for side necks	29/32
Type piston	Three-neck round-bottom flask

