



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz nitrogen round bottom flask (Schlenk flask) with NS stopcock, 100 ml, NS 14/23, NS stopcock 2.5 mm bore

€56.85
plus VAT &
Shipping

Product Images



Description

Nitrogen round bottom flasks made of DURAN® borosilicate glass 3.3 in various volumes. With NS conical ground joint and NS stopcock. With olive for tubing with an inner diameter of 10 mm. Suitable for the Schlenk technique.

Additional Information

No.	LZ-0730
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2730
EAN	4058072304706
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Pistons properties	with side necks with standard ground joint with NS stopcock/valve
Number of side necks	1
Neck size	not specified
Alignment side necks	Oblique 20°
Type cock / valve	NS cock with glass chick
Standard ground piston	14/23
Type piston	Flask round

