



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz nitrogen round bottom flask (Schlenk flask) with NS stopcock, 500 ml, NS 29/32, NS stopcock 2.5 mm bore

**€78.51**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Nitrogen round bottom flasks made of DURAN® borosilicate glass 3.3 in various volumes. With NS conical ground joint and NS stopcock. With olive for tubing with an inner diameter of 10 mm. Suitable for the Schlenk technique.

## Additional Information

---

No.	LZ-0735
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2735
EAN	4058072304751
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
Pistons properties	with side necks with standard ground joint with NS stopcock/valve
Number of side necks	1
Neck size	not specified
Alignment side necks	Oblique 20°
Type cock / valve	NS cock with glass chick
Standard ground piston	29/32
Type piston	Flask round

