



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Lenz nitrogen round bottom flask (Schlenk flask) with PTFE valve stopcock, 250 ml, NS 14/23, valve bore 2.5 mm

**€84.71**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

With NS conical ground joint and PTFE valve stopcock, with labeling field. Material: DURAN®

## Additional Information

---

No.	LZ-0740
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2740
EAN	4058072304805
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	250 ml
Pistons properties	with side necks with standard ground joint with NS stopcock/valve
Number of side necks	1
Neck size	not specified
Alignment side necks	Parallel
Type cock / valve	Valve with PTFE piston
Standard ground piston	14/23
Type piston	Flask round

