



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz nitrogen round-bottom flask (Schlenk flask) with PTFE valve stopcock, 50 ml, NS 29/32, valve bore 2.5 mm

€66.40
plus VAT &
Shipping

Product Images



Description

With NS conical ground joint and PTFE valve stopcock, with labeling field. Material: DURAN®

Additional Information

No.	LZ-0741
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2741
EAN	4058072304812
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
Pistons properties	with side necks with standard ground joint with NS stopcock/valve
Number of side necks	1
Neck size	not specified
Alignment side necks	Parallel
Type cock / valve	Valve with PTFE piston
Standard ground piston	29/32
Type piston	Flask round

