



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Lenz Nitrogen Template (Schlenk Tube) with NS stopcock, 25 ml, sleeve NS 14/23, NS stopcock 2.5 mm hole

**€57.75**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Schlenk tubes/nitrogen templates made of DURAN® borosilicate glass 3.3 in various volumes. With NS conical ground joint and NS stopcock. With valve bore: 2.5 mm With olive for tubes with an inner diameter of 10 mm. Suitable for Schlenk technique.

## Additional Information

---

No.	LZ-0747
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2747
EAN	4058072304874
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	25 ml
Pistons properties	with standard ground joint with NS stopcock/valve
Standard ground sleeve	14/23
Type cock / valve	NS cock with glass chick
Standard ground piston	14/23
Type piston	Swivel tube

