




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz nitrogen template (Schlenk tube) with NS
stopcock, 10 ml, core NS 14/23, NS stopcock 2.5
mm hole**

**€67.27
plus VAT &
Shipping**



Description

Schlenk tubes/nitrogen templates made of DURAN® borosilicate glass 3.3 in various volumes. With NS conical ground joint and NS stopcock. With valve bore: 2.5 mm With olive for tubes with an inner diameter of 10 mm. Suitable for Schlenk technique.

Additional Information

No.	LZ-0751
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2751
EAN	4058072304911
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	10 ml
Pistons properties	with standard ground joint with NS stopcock/valve
Standard ground core	14/23
Type cock / valve	NS cock with glass chick
Standard ground piston	14/23
Type piston	Swivel tube

