



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz nitrogen receiver (Schlenk tube) with PTFE  
valve stopcock, 25 ml, sleeve NS 14/23, valve  
bore 2.5 mm**

**€69.49  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Nitrogen receivers with NS conical ground joint and PTFE valve stopcock; made of DURAN® tube

## Additional Information

---

No.	LZ-0757
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2757
EAN	4058072304973
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	25 ml
Pistons properties	with standard ground joint with NS stopcock/valve
Standard ground sleeve	14/23
Type cock / valve	Valve with PTFE piston
Standard ground piston	14/23
Type piston	Swivel tube

