



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz two-neck pointed flask, 50 ml, MH NS 14/23, SH NS 14/23, DURAN®

€34.20
plus VAT &
Shipping

Product Images



Description

Two-neck pointed flask

According to DIN 12383 and in addition to the DIN standard, with slanted side neck (angle: 25°).
With standard ground joint, with inscription field.

Additional Information

No.	LZ-0622
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2622
EAN	4058072303624
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
Pistons properties	with side necks with standard grinding
Number of side necks	1
Standard ground sleeve	14/23
Alignment side necks	Oblique 25°
DM Piston	50 mm
Standard grinding for side necks	14/23
Type piston	Pointed flask

