



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz upright flask, narrow neck, capacity 4,000 ml, with beaded rim, DURAN®

**€112.17
plus VAT &
Shipping**

Product Images



Description

Flask with beaded rim Narrow neck design According to DIN EN ISO 4797 and in addition to the standard, with labeling field, Material: DURAN®

Additional Information

No.	LZ-2198
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4202
EAN	4058072319311
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	4000 ml
Neck size	narrow
DM Neck	50 mm
DM Piston	207 mm
Type piston	Plunger

