



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Lenz upright flask, wide neck, capacity 500 ml,  
with beaded rim, DURAN®**

**€152.00  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Flask with beaded rim Wide-necked version According to DIN EN ISO 4797 and in addition to the standard, with labeling field, Material: DURAN®

## Additional Information

---

No.	LZ-2205
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4209
EAN	4058072319380
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
Neck size	far
DM Neck	50 mm
DM Piston	103 mm
Type piston	Plunger

