



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**labtable laboratory hanging cabinet 1 glass
door, 1 shelf WxHxD: 600 mm x 740 mm x 350
mm**

**€300.80
plus VAT &
Shipping**

Product Images



Description

labtable laboratory hanging cabinets provide storage space for chemicals, laboratory materials and small equipment. For an optimal storage solution, we offer various designs of cabinets with glass and full doors, as well as open shelves.

- Made of melamine
- Doors equipped with dampers
- Made of laminated chipboard with a thickness of 18 mm
- Doors with ABS edge thickness of 2 mm

Additional Information

| | |
|-------------------------------|----------------------|
| No. | G-0069 |
| Manufacturer (Brand) | labtable |
| EAN | 4058072254872 |
| Transport temperature | Room temperature |
| Color | light gray |
| Material | Melamine |
| Wide | 600 mm |
| Height | 740 mm |
| Depth | 350 mm |
| Number of shelves/ grids | 1 |
| Number of doors | 1 |
| Page door suspension | left |
| Type Shipping / Delivery Cost | Special shipping fee |

