



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



labtable Laboratory Hanging Cabinet 2 glass doors, 1 shelf WxHxD: 900 mm x 740 mm x 350 mm

€438.30
plus VAT &
Shipping

Product Images



Description

labtable laboratory hanging cabinets provide storage space for chemicals, laboratory materials and small equipment. For an optimal storage solution, we offer various designs of cabinets with glass and full doors, as well as open shelves.

- Made of melamine
- Doors equipped with dampers
- Made of laminated chipboard with a thickness of 18 mm
- Doors with ABS edge thickness of 2 mm

Additional Information

No.	G-0072
Manufacturer (Brand)	labtable
EAN	4058072254902
Transport temperature	Room temperature
Color	light gray
Material	Melamine
Wide	900 mm
Height	740 mm
Depth	350 mm
Number of shelves/ grids	1
Number of doors	2
Page door suspension	left right
Type Shipping / Delivery Cost	Special shipping fee

