



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Drawer insert for trolley 60 cm, variable, transparent

€73.40
plus VAT &
Shipping

Product Images



Description

The drawer insert with variable division (3 lanes, 6 dividers) adapts to your needs.

Transparent plastic, dishwasher safe.

Dimensions: 500 x 385 x 55 mm (W x D x H)

Accessories for: Trolley with drawers (Order No. 2-7022 and 2-7023)

Undercarriage, 52 cm high (Order No. 2-7030 and 2-7031)

Additional Information

No.	2-7024
Manufacturer (Brand)	neoLab
EAN	4058072008505
Transport temperature	Room temperature
Color	White
Material	Polycarbonate (PC)
Wide	385 mm
Height	55 mm
Length	500 mm
Depth	500 mm
Number of compartments	3

