



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



RSG stainless steel wire mesh, 155 x 155 mm

€10.93 As low as
€8.41
plus VAT &
Shipping

Product Images



Description

Wire nets made entirely of heat-resistant stainless steel for burners, as a support for tribrachs or the like. The wire mesh is provided with a sturdy stainless steel sheet frame.
The particularly fine mesh center effectively prevents the flames from passing through.

Additional Information

No.	RG-1087
Manufacturer (Brand)	RSG
old neoLab article no.	1-6777
EAN	4058072338633
Transport temperature	Room temperature
Material	Stainless steel
Wide	155 mm
Length	155 mm
Protective shields and plates Features	with frame
Thickness	2.5 mm
Dimensions	155 x 155 mm

