



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Parafilm dispenser made of plexiglass, blue

€127.60
plus VAT &
Shipping

Product Images



Description

Plexiglas dispenser for parafilm standard rolls with 10 cm width.

Withhinged lid and cutter.

Film is stored neatly, cleanly and clearly.

With improved feeding of the parafilm and easily movable cutting knife.

Dimensions: 12 x 16 x 17 cm (W x L x H) with lid.

Additional Information

No.	3-1004
Manufacturer (Brand)	neoLab
EAN	4058072011864
Transport temperature	Room temperature
Color	transparent (blue)
Material	Polymethylmetacrylate (PMMA. Plexiglas)
Wide	120 mm
Form	angular
Height	170 mm
Length	160 mm
Depth	160 mm
Dispenser properties	with closure/lid
for width tape/ roll	up to 100 mm

