



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Storage box for earmuffs, transparent, midi

€45.60
plus VAT &
Shipping

Product Images



Description

Safety must be at hand! This is possible with the storage box, ideal for all common earmuffs.

Opening handle can be used for nameplate.

Labeled with international command symbol according to VGB 125.

Dust and splash-proof.

Easy to attach by sticking or screwing on.

Fastening material and labels included.

Quick removal.

Additional Information

No.	2-6060
Manufacturer (Brand)	neoLab
EAN	4058072069988
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky) black
Material	Acrylonitrile butadiene styrene (ABS)
Wide	125 mm
Height	225 mm
Length	236 mm
Depth	236 mm
suitable for	Earplugs/ Ear muffs
Dispenser properties	with closure/lid suitable for wall mounting
Capacity	variable (depending on vessel volume)

