



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Storage box for earmuffs, blue, midi

**€38.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Safety must be at hand! This is possible with the storage box, ideal for all common earmuffs.**

Opening handle can be used for name plate.

Labeled with international command symbol according to VGB 125.

Dust and splash-proof.

Easy to attach by sticking or screwing on.

Fastening material and labels included.

Quick removal.

## Additional Information

---

No.	2-6061
Manufacturer (Brand)	neoLab
EAN	4058072069995
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue black
Material	Acrylonitrile butadiene styrene (ABS)
Wide	125 mm
Height	225 mm
Length	236 mm
Depth	236 mm
suitable for	Earplugs/ Ear muffs
Dispenser properties	with closure/lid suitable for wall mounting

