



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Storage box f. Hair caps, transparent, midi

**€49.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Hair protection hoods always ready to hand!

Attach storage box in entrance areas to pharmaceutical production rooms, laboratories or clean rooms.

Opening handle can be used as a label.

Labeling with international mandatory signs according to VGB 125.

Dust and splash-proof.

Easy attachment , ABS material, functional design.

## Additional Information

---

No.	2-6063
Manufacturer (Brand)	neoLab
EAN	4058072004910
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky) black
Material	Acrylonitrile butadiene styrene (ABS)
Wide	125 mm
Height	225 mm
Length	236 mm
Depth	236 mm
suitable for	Hair protection hoods
Dispenser properties	with closure/lid suitable for wall mounting

