



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Storage box f. Hair caps, blue, midi

€34.30
plus VAT &
Shipping

Product Images



Description

Hair protection hoods always at hand!

Attach storage box in entrance areas to pharmaceutical production rooms, laboratories or clean rooms.

Opening handle can be used as a label.

Labeling with international mandatory signs according to VGB 125.

Dust and splash-proof.

Easy attachment , ABS material, functional design.

Additional Information

| | |
|-----------------------|---|
| No. | 2-6064 |
| Manufacturer (Brand) | neolab |
| EAN | 4058072004927 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | blue black |
| Material | Acrylonitrile butadiene styrene (ABS) |
| Wide | 125 mm |
| Height | 225 mm |
| Length | 236 mm |
| Depth | 236 mm |
| suitable for | Hair protection hoods |
| Dispenser properties | with closure/lid suitable for wall mounting |

