



neoLab® Glove dispenser orange

Product Images



Description

Orange glove dispenser for disposable glove boxes with 100 pieces content.

Protects the boxes from splashes.

Boxes are inserted from the side and the gloves are removed from the top. For quick identification of different glove sizes or materials.

Dispenser suitable for wall mounting.

Inner dimensions: 248 x 130 x 80 mm (L x D x H)

Additional Information

No.	1-6800
Manufacturer (Brand)	neoLabLine
EAN	4058072044404
Color	transparent (orange)
Material	Polymethylmetacrylat (PMMA) (Acrylglas)
Wide	90 mm
Height	140 mm
Length	250 mm
Depth	250 mm
Dispenser properties	für die Wandmontage geeignet
Maximum width box	130 mm
Maximum height box	80 mm
Maximum length box	250 mm

