



## neoLab® glove dispenser red

### Product Images

---



### Description

---

Red glove dispenser for disposable glove boxes with 100 pieces content.

Protects the boxes from splashes.

Boxes are inserted from the side and the gloves are removed from the top. For quick identification of different glove sizes or materials.

Dispenser suitable for wall mounting.

**Inner dimensions:** 248 x 130 x 80 mm (L x D x H)

## Additional Information

|                      |                                         |
|----------------------|-----------------------------------------|
| No.                  | 1-6802                                  |
| Manufacturer (Brand) | neoLabLine                              |
| EAN                  | 4058072044428                           |
| Color                | transparent (rot)                       |
| Material             | Polymethylmetacrylat (PMMA) (Acrylglas) |
| Wide                 | 90 mm                                   |
| Height               | 140 mm                                  |
| Length               | 250 mm                                  |
| Depth                | 250 mm                                  |
| Dispenser properties | für die Wandmontage geeignet            |
| Maximum width box    | 130 mm                                  |
| Maximum height box   | 80 mm                                   |
| Maximum length box   | 250 mm                                  |

