



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Universal Stainless Steel Glove Dispenser "Multismart"

**€82.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

neoLab Universal **glove dispenser** "Multismart" for 3 boxes

This dispenser made of high-quality stainless steel is easily adjustable for both the depth of the packaging and the width of the box. The height of the packages is even variable without any changes. The dispenser is manufactured in stainless steel due to the hygiene issue, other designs on request! Made of 2 parts, with closed back wall. Intended for wall mounting as standard, mounting on equipment rail and round bars is also possible as an option. Screws and dowels included.

Width variants: 200 - 260 mm

Depth variants: 40 - 95 mm

Total height: 39 cm

## Additional Information

---

No.	1-0680
Manufacturer (Brand)	neoLab
EAN	4058072115845
Color	silver
Material	Stainless steel
Wide	120 mm
Height	390 mm
Length	250 mm
Depth	250 mm
Dispenser properties	suitable for wall mounting
Maximum width box	130 mm
Maximum height box	90 mm
Maximum length box	230 mm

