



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® High quality glove dispenser stainless steel, for 1 pack

€42.10
plus VAT &
Shipping

Product Images



Description

With the simple elegance of stainless steel!

Glove dispenser for almost all glove or cosmetic wipes made of high quality stainless steel, matt brushed surface. Quick to remove from the wall and reattach without tools thanks to the keyhole system.

Order No. 1-6804

Dispenser for 1 box

Suitable for all dispenser boxes with maximum dimensions 250 x 132 x 90 mm (L x W x H)

Additional Information

No.	1-6804
Manufacturer (Brand)	neoLabLine
EAN	4058072044442
Color	silver
Material	Stainless steel
Wide	135 mm
Height	95 mm
Length	250 mm
Depth	250 mm
Dispenser properties	suitable for wall mounting
Maximum width box	132 mm
Maximum height box	90 mm
Maximum length box	250 mm

