



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® High quality glove dispenser made of stainless steel, flexible for 1 box with max. 250x132x60-90 mm

**€53.30
plus VAT &
Shipping**

Product Images



Description

Dispenser for 1 carton

flexible: in the bottom of the dispenser there is a tab that can be bent by the user, so cartons with different heights can be inserted and always securely fastened.

suitable for all dispenser cartons with maximum dimensions 250 x 132 x 60 to 90 mm (L x W x H).

Additional Information

No.	1-6806
Manufacturer (Brand)	neoLabLine
EAN	4058072092139
Color	silver
Material	Stainless steel
Wide	135 mm
Height	75 mm
Length	250 mm
Depth	250 mm
Dispenser properties	suitable for wall mounting
Maximum width box	132 mm
Maximum height box	75 mm
Maximum length box	250 mm

