



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® glove dispenser made of stainless steel for 1 box with max. 250x135x95 mm

**€17.50
plus VAT &
Shipping**

Product Images



Description

Dispenser for 1 carton. Suitable for all dispenser boxes with maximum dimensions 250 x 132 x 90 mm (L x W x H).

Additional Information

No.	1-6808
Manufacturer (Brand)	neoLabLine
EAN	4058072092115
Color	silver
Material	Stainless steel
Wide	135 mm
Height	95 mm
Length	250 mm
Depth	250 mm
Dispenser properties	suitable for wall mounting
Maximum width box	135 mm
Maximum height box	95 mm
Maximum length box	250 mm

