



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® glove dispenser made of stainless  
steel for 3 boxes with max. 258x392x96 mm**

**€35.00  
plus VAT &  
Shipping**

# Product Images



## Description

---

Order No. 1-6809. Dispenser for 3 cartons. Suitable for all dispenser cartons with maximum dimensions 250 x 130 x 90 mm (L x W x H). Punched marking for 1 carton each of size L, M and S.

## Additional Information

---

No.	1-6809
Manufacturer (Brand)	neoLabLine
EAN	4058072092122
Color	silver
Material	Stainless steel
Wide	258 mm
Height	95 mm
Length	392 mm
Depth	392 mm
Dispenser properties	suitable for wall mounting
Capacity	3
Maximum width box	135 mm
Maximum height box	95 mm
Maximum length box	250 mm

