



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Petri dish rack, stackable

**€24.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Stackable Petri dish trays for incubators.  
Holes in the bottom allow optimal circulation of air.  
Can also be used as a stacking rack on the lab bench.  
Trays are made of white polypropylene with rim.  
Each tray has 2 labeling fields for sample identification.

Loading: 5 x 100 mm Ø or 15 x 60 mm Ø

**Dimensions:** 240 x 250 x 30 mm (W x D x H)

## Additional Information

---

No.	2-1306
Manufacturer (Brand)	neoLab
EAN	4058072056001
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polycarbonate (PC)
Wide	251 mm
Height	31 mm
Length	242 mm
Depth	242 mm
suitable for	Petri dishes
Dispenser properties	stackable / interlockable
Type frame	Petri dish rack
for vessels with DM up to	60.0 mm 100.0 mm

