



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Storage Box Maxi with Insert for 12 Safety Goggles

**€127.10**  
plus VAT &  
Shipping

### Product Images



## Description

---

The ideal glasses dispenser with inserts for 12 safety glasses. Space-saving and clearly arranged, the glasses are stored.

Quick access to the glasses due to the grid insert.

For production entrances, laboratories or other public areas.

Affixed mandatory sign according to

BGV A 8.

Scope of delivery: fixing material for wall mounting

## Additional Information

---

No.	1-1422
Manufacturer (Brand)	neoLab
EAN	4058072032890
autoclavable at 121°C	No
Color	blue black
Material	Acrylonitrile butadiene styrene (ABS)
Wide	200 mm
Height	315 mm
Length	236 mm
Depth	236 mm
suitable for	Protective goggles
Dispenser properties	with closure/lid suitable for wall mounting
Capacity	12

