



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Storage box for safety goggles, blue, midi

**€37.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Safe and clean storage box for all safety glasses.

Directly at the workplace for employees and also for visitors.

Quick and easy handling due to **quick-open/quick-close closure**.

Opening handle can be used for name tag. Labels and fastening material included.

Dust- and splash-proof.

Simple attachment by gluing or screwing on.

Safety due to the strong hygienic material ABS in signal color with international symbols according to VGB 125.

Quick removal, functional design.

User-friendly, structured surface.

## Additional Information

---

No.	2-6056
Manufacturer (Brand)	neoLab
EAN	4058072004873
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue black
Material	Acrylonitrile butadiene styrene (ABS)
Wide	125 mm
Height	236 mm
Length	225 mm
Depth	225 mm
suitable for	Protective goggles mouthguard gloves
Dispenser properties	with closure/lid suitable for wall mounting

