



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Storage box f. Gloves, blue, midi

€39.30
plus VAT &
Shipping

Product Images



Description

Safe and clean storage box for all gloves.

Directly at the workplace for employees and also for visitors.

Quick and easy handling due to **quick-open/quick-close closure**.

Opening handle can be used for name tag. Labels and fastening material included.

Dust- and splash-proof.

Simple attachment by gluing or screwing on.

Safety due to the strong hygienic material ABS in signal color with international symbols according to VGB 125.

Quick removal, functional design.

User-friendly, structured surface.

Additional Information

No.	2-6058
Manufacturer (Brand)	neoLab
EAN	4058072004897
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue black
Material	Acrylonitrile butadiene styrene (ABS)
Wide	125 mm
Height	236 mm
Length	225 mm
Depth	225 mm
suitable for	Protective goggles mouthguard gloves
Dispenser properties	with closure/lid suitable for wall mounting
Capacity	variable (depending on vessel volume)

