



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® safety glasses dispenser acrylic for up to 10 glasses, 323 x 178 x 98 mm

€122.00
plus VAT &
Shipping

Product Images



Description

Ideal for storing visitor goggles!

Up to 5 or 10 pairs of protective goggles are stored in this brown acrylic container, quickly at hand and protected from dust.

The goggles are inserted from above and removed through the slot at the bottom.

Mounting material for wall mounting is included.

Goggles are not included.

Additional Information

No.	2-4166
Manufacturer (Brand)	neoLabLine
EAN	4058072002336
Transport temperature	Room temperature
Color	transparent (brown)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	98 mm
Height	323 mm
Length	178 mm
Depth	178 mm
suitable for	Goggles / safety glasses
Dispenser properties	mobile with closure/lid suitable for wall mounting

