



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® safety glasses dispenser acrylic for up to 5 glasses, 198 x 178 x 98 mm

**€118.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Ideal for storing visitor goggles!

Up to 5 or 10 pairs of protective goggles are stored in this brown acrylic container, quickly at hand and protected from dust.

The goggles are inserted from above and removed through the slot at the bottom.

Mounting material for wall mounting is included.

Goggles are not included.

## Additional Information

---

No.	2-4167
Manufacturer (Brand)	neoLabLine
EAN	4058072002343
Transport temperature	Room temperature
Color	transparent (brown)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	98 mm
Height	198 mm
Length	178 mm
Depth	178 mm
suitable for	Goggles / safety glasses
Dispenser properties	mobile with closure/lid suitable for wall mounting

