



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Dispenser for liquid soap and hand sanitiser with push button, for wall mounting, 1000 ml

€23.93
plus VAT &
Shipping

Product Images



Description

neoLab dispenser for hand disinfectant with push button, for wall mounting, 1000 ml

Thanks to its flexibility, the neoLab wall-mounted dispenser can be easily and quickly attached to the desired position. Fastening material enclosed. To refill the disinfectant, the integrated tank can be removed, making refilling **child's** play. The tank holds a volume of 1000ml.

The desired amount of hand disinfectant (1.8 ml per dosing operation) is dispensed precisely via a pressure switch. **No batteries** are necessary for use!

Additional Information

No.	1-0697
neoGreen	No
Manufacturer (Brand)	neoLab
EAN	4058072122874
autoclavable at 121°C	No
Color	white black
Material	Plastic
sterile	No
Form	rectangular
Height	260 mm
Depth	132 mm
suitable for	Liquid soap hand disinfectant
Energy supply	Battery power supply
Type operation	Touchless

