



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Dispenser for hand disinfectant with infrared sensor, 350 ml

Special Price

**€17.50** was

**€19.09**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

### neoLab Dispenser for hand disinfectant with infrared sensor, 350 ml.

The neoLab dispenser can be easily used at any location due to its mobility. The contactless hand disinfectant dispensing is precise. The reservoir has a volume of up to 350ml and can be refilled easily. The dispensing volume is preset and cannot be adjusted individually.

**Note: Batteries are not included!**

## Additional Information

---

No.	1-0698
Manufacturer (Brand)	neoLab
EAN	4058072122881
diameter	65 mm
Color	whiteltransparent (milky)
Material	Plastic
Form	cylindrical
Height	230 mm
suitable for	Liquid soap hand disinfectant
Energy supply	Battery
Type operation	Touchless

