



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Hygomat classic overshoe dispenser, automatic

€447.00
plus VAT &
Shipping

Product Images



Description

The **overshoe dispenser** HYGOMAT Classic is an excellent solution for all who **want to save time during the annoying changing**. With the HYGOMAT, putting on overshoes is 4x faster! A time saving that really pays off. At the same time, this overshoe automat effectively prevents cross-contamination and accidents. Because overshoes ensure hygiene and a safe step in hygiene-sensitive areas.

- can be set up anywhere and very easily
- no power supply necessary
- very user-friendly, no bending, easy and fast handling
- four times faster than manual overshoe donning
- with viewing window, so that the filling level can be seen at any time
- refilled with overshoes in a few seconds
- Cleanroom certification ISO class 5* according to EN ISO 14644-1

Additional Information

No.	2-0080
Manufacturer (Brand)	neoLab
EAN	4058072161910
Color	White Grey
Material	Acrylonitrile Butadiene Styrene (ABS) Stainless Steel
Wide	315 mm
Height	230 mm
Length	695 mm
Depth	695 mm
suitable for	Overshoes
Dispenser properties	mobile

