



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Hygomat comfort overshoe dispenser, automatic

€536.10
plus VAT &
Shipping

Product Images



Description

The **overshoe dispenser** HYGOMAT COMFORT is an excellent solution for all those who **want to save time during the tiresome changing process**. With the HYGOMAT, putting on overshoes is 4x faster! A time saving, which is really worthwhile. At the same time, this overshoe automat effectively prevents cross-contamination and accidents. Because overshoes ensure hygiene and a safe step in hygiene-sensitive areas.

- can be set up anywhere and very easily
- no power supply necessary
- very user-friendly, no bending, easy and fast handling
- **with handle to hold - for even more comfortable donning**
- four times faster than manual overshoe donning
- with viewing window, so that the filling level can be seen at any time
- refilled with overshoes in a few seconds
- Cleanroom certification ISO class 5* according to EN ISO 14644-1

Additional Information

No.	2-0081
Manufacturer (Brand)	neoLab
EAN	4058072155254
Color	White Grey
Material	Acrylonitrile Butadiene Styrene (ABS) Stainless Steel
Wide	300 mm
Height	730 mm
Length	740 mm
Depth	740 mm
suitable for	Overshoes
Dispenser properties	mobile

