



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Moonlab® Small parts dispenser small, acrylic, yellow, (LxWxH) 200x140x190 mm

€43.90
plus VAT &
Shipping

Product Images



Description

The Moonlab small parts dispensers help you to organize various small parts clearly and in a space-saving manner.

Filling from above via removable lid.

Removal from below, via an inclined chute with hinged lid.

There are 2 suspension eyes on the back for wall mounting.

Pipette tips, Pasteur pipettes or centrifuge tubes, for example, can be stored in the dispensers.

The dispensers are made of acrylic plastic and are therefore extremely durable.

Color: yellow

Size: small

Additional Information

No.	4-0079
Manufacturer (Brand)	Moonlab
EAN	4058072086510
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	yellow
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	140 mm
Height	190 mm
Length	200 mm
Depth	200 mm
Dispenser properties	with closure/lid suitable for wall mounting
Number of compartments	1

