



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Moonlab® Small parts dispenser large, acrylic,
black, (LxWxH) 200x140x275 mm, 1 pc**

**€51.00
plus VAT &
Shipping**

Product Images



Description

The Moonlab small parts dispensers help you to organize various small parts clearly and space-saving.

Filling from the top through removable lid.

Removal from the bottom, via a slanted chute with hinged lid.

For wall mounting, there are 2 hanging eyes on the back.

For example, pipette tips, Pasteur pipettes or centrifuge tubes can be stored in the dispensers.

The dispensers are made of acrylic plastic and are therefore extremely durable.

Color: black

Size: large

Additional Information

No.	4-0083
Manufacturer (Brand)	Moonlab
EAN	4058072087784
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	140 mm
Height	275 mm
Length	200 mm
Depth	200 mm
Dispenser properties	with closure/lid suitable for wall mounting
Number of compartments	1

