



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® assortment box (PS), 20 compartments

**€18.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Assortment box with hinged lid made of crystal-clear polystyrene, for clear storage of small parts in the laboratory.

2 rows with 10 openings each,  
each 77 x 27 x 38 mm deep.

Temperature resistant from -40°C to +70°C

Outer dimensions: 300 x 170 x 45 mm (L x W x H)

## Additional Information

---

No.	2-1540
Manufacturer (Brand)	neoLab
EAN	4058072057305
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polystyrene (PS)
Wide	170 mm
Height	45 mm
Length	300 mm
Depth	300 mm
suitable for	variable applications
Dispenser properties	stackable/assemblable with closure/lid
Number of compartments	20

