



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Acrylic dispenser with 3 compartments, 14 x 16 x 21.5 cm

€217.60
plus VAT &
Shipping

Product Images



Description

The ideal dispenser and storage box for a variety of small parts, such as pipette tips, disposable syringes, reaction vessels, etc., metal or glass for laboratory and service.

The box is made of 3 mm thick, durable acrylic plastic.

Loading is via the lid, removal via the front openings.

In addition to easy accessibility, the parts are clearly stored and protected against dust.

Additional Information

No.	2-1516
Manufacturer (Brand)	neoLabLine
EAN	4058072057206
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	210 mm
Height	215 mm
Length	160 mm
Depth	160 mm
suitable for	variable applications
Dispenser properties	with closure/lid
Number of compartments	3

