



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Acrylic dispenser with 4 compartments, 29 x 16 x 21.5 cm

**€254.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The ideal dispenser and storage box for a variety of small parts, such as pipette tips, disposable syringes, reaction vessels, etc., metal or glass for laboratory and service.

The box is made of 3 mm thick, durable acrylic plastic.

Loading is via the lid, removal via the front openings.

In addition to easy accessibility, the parts are clearly stored and protected against dust.

## Additional Information

---

No.	2-1517
Manufacturer (Brand)	neoLabLine
EAN	4058072057213
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	290 mm
Height	215 mm
Length	160 mm
Depth	160 mm
suitable for	variable applications
Dispenser properties	with closure/lid
Number of compartments	4

