



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Assortment Box Standard Shape, Transparent with Transparent Lid

**€12.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Practical neoLab system box made of ABS plastic with attached hinged lid and sturdy snap closure. The boxes are stackable so that slipping can be prevented. – available in various sizes and colors –

**Suitable for two inserts of your choice**

## Additional Information

---

No.	2-0051
Manufacturer (Brand)	neoLab
EAN	4058072089443
autoclavable at 121°C	No
Color	whiteltransparent (milky)
Material	Acrylonitrile butadiene styrene (ABS)
Wide	126 mm
Height	37 mm
Length	190 mm
Depth	190 mm
Dispenser properties	stackable/collapsible with removable insert with closure/lid
Number of compartments	1

