



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® assortment box high form, white with transparent lid

€13.60
plus VAT &
Shipping

Product Images



Description

Practical neoLab system box made of ABS plastic with attached hinged lid and sturdy snap closure. The boxes are stackable so that slipping can be prevented. - available in various sizes and colors -

Suitable for three inserts of your choice

Additional Information

No.	2-0052
Manufacturer (Brand)	neoLab
EAN	4058072089450
autoclavable at 121°C	No
Color	transparent (milky)
Material	Acrylonitrile butadiene styrene (ABS)
Wide	126 mm
Height	96 mm
Length	190 mm
Depth	190 mm
Dispenser properties	stackable/collapsible with removable insert with closure/lid
Number of compartments	1

