



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® System box transp. PS, 230 x 150 x 130 mm

€13.30 As low as
€10.50
plus VAT &
Shipping

Product Images



Description

System boxes with proven features:

Crystal clear polystyrene gives immediate overview of contents & quantity.

Large viewing opening, high stackability.

Wide label field, excellent stability.

Accessories: Labels for System Box (Order No. 2-1295)

Additional Information

No.	2-1291
Manufacturer (Brand)	neoLab
EAN	4058072055912
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polystyrene (PS)
Wide	230 mm
Height	150 mm
Length	130 mm
Depth	130 mm
suitable for	variable applications
Dispenser properties	stackable/pluggable with labeling field
Number of compartments	1

