



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® laboratory stool with castors and foot
ring, PU comfort foam seat, height-adjustable
59-83 cm**

**€401.40
plus VAT &
Shipping**

Product Images



Description

Laboratory stool with five-legged base in chromed steel with casters or oscillating glides. The seat is made of comfort polyurethane foam, easy to decontaminate, suitable for clean rooms class S1, S2.

Chromed column, infinitely variable seat height adjustment by means of gas spring (depth suspension) from 55 to 83 cm.

Aluminum foot ring height adjustable.

GS tested.

Additional Information

No.	1-1074
Manufacturer (Brand)	neoLab
EAN	4058072031800
Color	silver
Material	Steel chrome plated
Chairs, stools and standing aids features	suitable for cleanroom with footrest/ring
Figure	2
Height adjustment	59 to 83 cm
Material frame	Steel chrome plated
Type foot	with rollers
Frame color	silver
Color seat / backrest	black
Material rollers	Polypropylene (PP)
Material seat / backrest	Polyurethane (PU)
Type chair	Swivel stool

