



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® cleanroom chair PU foam black,
height adj. 52-65 cm,with castors**

**€392.70
plus VAT &
Shipping**



Description

Comfortable chairs specially designed for the cleanroom, durable and functional.

Seat and back in soft black comfort polyurethane foam (4 cm thickness), easy to decontaminate. The soft polyurethane foam is resistant to cuts, flying sparks, light alkalis, acids and lubricants.

- Seat height adjustment via gas spring
- depth adjustable
- chromed five-legged base with castors and glides
- chrome-plated pedestal, gas spring height adjustment
- Depth suspension

GS tested

Additional Information

| | |
|---|------------------------|
| No. | 1-1070 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072031763 |
| Color | silver |
| Material | Steel chrome plated |
| Chairs, stools and standing aids features | suitable for cleanroom |
| Figure | 1 |
| Frame color | silver |
| Color seat / backrest | black |
| Height adjustment | 52 to 65 cm |
| Material frame | Steel chrome plated |
| Material rollers | Polypropylene (PP) |
| Material seat / backrest | Polyurethane (PU) |
| Type foot | with rollers |
| Type chair | Swivel / office chair |

