



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## neolab® laboratory chair with castors, PU seat black, 45-58 cm

**€434.40**  
**plus VAT &**  
**Shipping**



## Description

---

### Comfortable chairs specially designed for the cleanroom, durable and functional.

Seat and back in soft blue comfort polyurethane foam (4 cm thickness), easy to decontaminate. The soft polyurethane foam is resistant to cuts, flying sparks, light alkalis, acids and lubricants.

- Seat height adjustment via gas spring
- depth adjustable
- chromed five-legged base with castors and glides
- chrome-plated pedestal, gas spring height adjustment
- Depth suspension

GS tested

## Additional Information

---

No.	1-1142
Manufacturer (Brand)	neoLab
EAN	4058072032050
Color	silver
Material	Steel chrome plated
Chairs, stools and standing aids features	with adjustable seat inclination
Figure	1
Height adjustment	48 to 58 cm
Material frame	Steel chrome plated
Type foot	with rollers
Frame color	silver
Color seat / backrest	black
Material rollers	Polypropylene (PP)
Material seat / backrest	Polyurethane (PU)
Type chair	Swivel / office chair

